



# VAP 10 Anaesthesia Vaporizer



**Halothane**



**Isoflurane**



**Sevoflurane**

**Vap 10 Vaporizer is designed to deliver accurate concentration under varying conditions of flow rate (even low flow) and temperatures**

- Available for Halothane, Isoflurane, Sevoflurane
- Supports Quick filling, Pour filling, Key filling
- VAP 10 Vaporizer compatible with Selectatec<sup>®</sup> with interlock, Cagemount, Drager Plug In, Drager Auto Exclusion mount
- The flexible design of VAP 10 Vaporizer compatible with other Manufacturers Anaesthetic workstations.
- Low Body Weight
- Maintenance Free

## Technical Specification :

DIMENSIONS	HEIGHT	WIDTH	DEPTH	APPROXIMATE WIEGHT
Selectatec® with interlock	240mm	120mm	200mm	5.2kg
Cagemount	227mm	132mm	187mm	5.0kg
Drager Plug In	227mm	132mm	204mm	5.2kg



### Output Volume accuracy

Accuracy  $\pm$  20% of scale Readings (or)  $\pm$  2% of Volume % (Whichever is greater)

### Capacity

Agent Capacity - 300 ml (with dry wick)

### Temperature Range

Operating Temperature Range - 15°C to 35°C

### Flow Range

Operating Flow range - 0.2 lpm to 15 lpm

### Pressure Range

Operating Pressure Range - 0 kpa to 5 kpa



### Safety Features


- Vaporizers are colour coded
- Clear agent level Indicator
- Actuating only one Vaporizer at a time
- MRI Compatible
- Tamper Proof Labels

### Vaporizer Type

	Sevoflurane 8% Filler Type			Isoflurane 5% Filler Type			Halothane 5% Filler Type		
	Pour fill	Key fill	Quick fill	Pour fill	Key fill	Easy fill	Pour fill	Key fill	Easy fill
Selectatec® Compatible with interlock	✓	✓	✓	✓	✓	✓	✓	✓	✓
Cagemount Taper 23 mm	✓	✓	✓	✓	✓	✓	✓	✓	✓
Drager Plug- In	✓	✓	✓	✓	✓	✓	✓	✓	✓
Drager Auto Exclusion	✓	✓	✓	✓	✓	✓	✓	✓	✓

### Typical performance graph

#### Variation of output



Quick fill drain adapter



Key fill adapter

### Backbar mounting System



Selectatec® with interlock



Cagemount Taper 23mm



Drager Plug In Compatible

